

Materials:

- Prepared slides
- Microscope
- Colored pencils

Procedure:

1. Pick up a prepared slide and write the name of the specimen in the results section below.
2. View the specimen under the microscope at low power.
3. Draw what you see under the microscope in the results section.
4. View the same specimen under medium power.
5. Draw what you see in the results section.
6. Repeat for high power and draw what you see.
7. Repeat steps 1-6 with a second prepared slide.




Results:

Name of Specimen _____

Low

Medium

High






Results:

Name of Specimen _____

Low

Medium


High



Results:


1. Name of item 1: _____
2. Draw what you observe using color.
3. Label the power with which you made your observations.
4. Describe your observation (what did you see)?

Power level _____



5. Name of item 2: _____
6. Draw what you observe using color.
7. Label the power with which you made your observations.
8. Describe your observation (what did you see)?

Power level _____



9. Name of item 3: _____
10. Draw what you observe using color.
11. Label the power with which you made your observations.
12. Describe your observation (what did you see)?

